

# ADDITIVE MANUFACTURING AND ASSEMBLY INSTRUCTIONS

## *ProGS-19*

Protective Glasses adjusted by sealing



## 1. Product Data

This is a protection device, called *ProGS-19*, totally developed for Open Source purpose. This device have been design on 20/03/2020 by Pablo Rodríguez González, Mechanical Design engineer, specialized on R&D Engineering and Design for manufacturing like Additive, Plastic extrusion, injection, metal machining parts and Composite Materials.

I invite you to contact me and help me with your feedback about this 1.0 version *ProGS-19* update through LinkedIn:

<https://www.linkedin.com/in/pablo-rodriguez-mechanical-engineer/>

This device has not any commercial purpose, it is only designed for helping the necessity of every medical worker during the complex context f the COVID-19 in Europe and world-wide. The ergonomic parameters have been followed to design for western complexion people

Due to this situation, the non-commercial use of this *ProGS-19* device is completely free for any MAKER technician or any person who wants to help the society around the world in this difficult situation.

You can download the file for free in the next link:

[\*ProGS-19 Download\*](#)

## 2. Bill of Material (BOM)

At first, we need a Bill of Materials to prepare in order to open the manufacturing process:

- 3D Printer



- PLA Filament



- Elastic strap (75cm at least).

Have a look to the link:

[https://www.amazon.es/Fellowes-Cristal-Portadas-encuadernaci3n-transparente/dp/B002BARBEK/ref=sr\\_1\\_1?\\_mk\\_es\\_ES=%C3%85M%C3%85%C5%BD%C3%95%C3%91&dchild=1&keywords=lamina+pvc+transparente&qid=1585321767&sr=8-1](https://www.amazon.es/Fellowes-Cristal-Portadas-encuadernaci3n-transparente/dp/B002BARBEK/ref=sr_1_1?_mk_es_ES=%C3%85M%C3%85%C5%BD%C3%95%C3%91&dchild=1&keywords=lamina+pvc+transparente&qid=1585321767&sr=8-1)

- Transparent colour laminar PVC (180-200 micrometer thickness).


Have a look to the link:

[https://www.amazon.es/LYTIVAGEN-El3stico-Fabricaci3n-Manualidades-Abalorios/dp/B07TF4GVYQ/ref=sr\\_1\\_1?\\_mk\\_es\\_ES=%C3%85M%C3%85%C5%BD%C3%95%C3%91&dchild=1&keywords=lytivagen+cordon+elastico&qid=1585321866&sr=8-1](https://www.amazon.es/LYTIVAGEN-El3stico-Fabricaci3n-Manualidades-Abalorios/dp/B07TF4GVYQ/ref=sr_1_1?_mk_es_ES=%C3%85M%C3%85%C5%BD%C3%95%C3%91&dchild=1&keywords=lytivagen+cordon+elastico&qid=1585321866&sr=8-1)






### 3. Manufacturing Instructions

Please, consider that the exact additive parameters will depend of the 3D printer used. In this case, we use both Creality-Ender3 and custommed Prusa i3Hephestos.

#### Step 1.



## Step 2.

 <b>Infill</b> <span>▼</span>	
<i>Infill Density</i>	<input type="text" value="12"/> %
Infill Line Distance	<input type="text" value="5.0"/> mm
<i>Infill Pattern</i>	<input type="text" value="Grid"/> <span>▼</span>
Infill Line Multiplier	<input type="text" value="1"/>
 <b>Material</b> <span>▼</span>	
Printing Temperature	<input type="text" value="200"/> °C
<i>Build Plate Temperature</i>	<input type="text" value="50"/> °C
Build Plate Temperature Initial Layer	<input type="text" value="50"/> °C
Flow	<input type="text" value="100"/> %
Enable Retraction	<input checked="" type="checkbox"/>
Retract at Layer Change	<input type="checkbox"/>
Retraction Distance	<input type="text" value="5"/> mm
Retraction Speed	<input type="text" value="45"/> mm/s
 <b>Speed</b> <span>ⓘ</span> <span>▼</span>	
<i>Print Speed</i>	<input type="text" value="60"/> mm/s
Enable Jerk Control	<input type="checkbox"/>
 <b>Travel</b> <span>&lt;</span>	
 <b>Cooling</b> <span>▼</span>	

### Step 3.

Profile

Standard Quality - 0.2mm



Q search settings



Retract at Layer Change

Retraction Distance

5 mm

Retraction Speed

45 mm/s



**Speed**



Print Speed



60 mm/s

Enable Jerk Control



**Travel**



**Cooling**



Enable Print Cooling

Fan Speed



100 %



**Support**



Generate Support



**Build Plate Adhesion**



Build Plate Adhesion Type



Skirt



Skirt Line Count



2

Skirt Distance



4 mm



**Dual Extrusion**



**Special Modes**




**Experimental**



Step 4.


Printing the file called: “ *Frame\_ProGS\_CV19\_v1\_1* ”




Step 5.

Printing the file called: “ *ClipWdow\_ProGS\_CV19* ”








Step 6.

Prepare all the components to develop the mechanical assembly.

Cut the Laminar PVC using the boundary shape of the file called "[Frame\\_ProGS\\_CV19\\_v1\\_1](#)", please look at the picture below:



Step 7.



Step 8.




#### 4. Complementary adjustment


Once you assemble the whole product, It will be interesting trying to improve the adjustment of the *ProGS-19* to the face complexion. It will be possible using an adhesive Isolation tape like it is showing in the following picture:



Then you will make a smoother surface finishing in the parts which are in physical contact with the person skin. Even you will improve the sealing of the *ProGS-19* to avoid any bacterial risk. See the picture below:



As you can see, the *ProGS-19* device design is provided of two little holes in order to combine the device with any standard glasses and use them simultaneously. Please have a look to the picture below:



Furthermore, the *ProGS-19* device can be sterilized manually using any Chloride or Alcoholic disinfectant gel. We use the following gel you can see in this picture:



ALCOHOLICA  
DE MANOS

©CARACTERÍSTICAS: Loción higienizante de manos sin agua, formulado con agentes bactericidas. Ideal para la limpieza e higiene de las manos en aquellas situaciones donde se recomienda la utilización de geles hidroalcohólicos por fricción.  
MODO DE EMPLEO: Dosificar una pequeña cantidad de producto en el hueco de la mano y frotar ambas manos hasta que se absorba y se evapore.  
©CARACTERÍSTICAS: Loção higienizante das mãos sem água formulado com agentes bactericidas. Ideal para a limpeza das mãos naquelas situações onde se recomenda a utilização dos geles hidroalcohólicos pela fricção.  
INSTRUÇÕES DE USO: Aplicar um pouco de produto na cavidade da mão e esfregando até que estiver absorvido e evaporado.  
INGREDIENTES: Alcohol denat, Aqua, Isopropyl alcohol, Glycerin, Hidroxipropyl Methylcellulose, Benzalkonium Chloride.



PELIGRO



## 5. Commentary

Finally, we necessary clarify the use of the *ProGS-19* is completely free purpose. Although, the commercial use of this technical research is totally forbidden.

Let's share this resource with everybody to help medical workers in hospitals and try to reduce the future impact of the Coronavirus.

Please, personally I ask you to collaborate with us manufacturing and sending every ProGS-19 protective glasses to the hospitals of your town.

Thank you in advance for your help and comprehension.

